nopa® Mobil

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MOBIL WITH ISOLATION TRANSFORMER**Description XU 800/41**

With the isolation transformer, the maximum device leakage current is reduced to a standard-conforming value.

When several electronic devices are used on one mobil, the individual devices are often supplied over a multiple socket strip, which is fixed to the mobil. The combined ground leakage current of these devices, then, can add up to a value, which is significantly higher than the allowed limit. If the protective conductor in the feed cable of the multiple socket strip is interrupted, a dangerous ground leakage current can flow through the patient to ground. When an isolation transformer is used, the ungrounded power supply of the individual devices prevents the ground leakage current from exceeding the allowed limit. This guarantees maximum safety for the patient, physicians and care personnel.

The isolation transformer can be installed under the bottom plate of the mobil XU 800/30 and XU 800/40.

Installation can only be carried out in the factory!

The isolation transformer is tested and certified according to IEC 601-1EN 60601-1, with a CE mark. It is conditionally short-circuit-proof, according to IEC 742 and IEC 601-1.

Technical data:

XU 800/41	Mobil with Isolation transformer
Power:	1600 VA
Protection class:	1
Inrush current limited to approx.:	20 A
Insulation voltage:	4 kV
Fuse protection:	primary 2-phase, secondary 1-phase

Including:

- central ON/OFF switch on the right E-column
- cold device socket strip (8 sockets), each with a POAG plug and fuse per socket, mounted on the left E-column
- set of rolls, Ø125 mm, conductive
- optionally with IsoMax insulation monitor XU 850/51 see catalogue page XU 850/51