

FLEXIBLE NASOPHARYNGOSCOPE



XQ 510/30

Das nopa Nasopharyngoskop wird komplett mit Aufbewahrungskoffer und Drucktester geliefert.
The nopa nasopharyngoscope will be delivered in a case and with a pressure testing device.

Beschreibung und technische Daten auf der Rückseite
Description and technical data on reverse side

**Flexibles Nasopharyngoskop****Beschreibung XQ 510/30**

Diese neue Generation des flexiblen Nasopharyngoskopes von nopa wurde unter Berücksichtigung der neuesten Produktionsmethoden und Technologien in enger Zusammenarbeit mit klinischen Anwendern entwickelt. Die innovativen Besonderheiten und die herausragende optische Qualität garantieren dem Benutzer ein optimales Bild.

Die neue Generation der Bildübertragungs-Technik liefert brillante Detailaufnahmen sowohl in den Ecken eines Bildes als auch in der Bildmitte. Das Erkennen kleinster Strukturen wird dadurch erheblich erleichtert. Die Wiedergabe der natürlichen Farben ermöglicht dem Benutzer eine bessere Diagnostik.

Das nopa Nasopharyngoskop bietet ein Höchstmaß an Hygiene.

Technische Daten:**XQ 510/30 nopa Nasopharyngoskop**

Bildwinkel:	85°
Tiefenschärfe:	2,5 - 5 mm
Außendurchmesser:	3,6 mm
Außendurchmesser distales Ende:	3,6 mm
Nutzlänge:	300 mm

Schwenkbereich des distalen Endes

Aufwärts:	125°
Abwärts:	125°

Flexible Nasopharyngoscope**Description XQ 510/30**

This new generation of the flexible nasopharyngoscope from nopa has been newly developed in close cooperation with clinical users, with an eye to production methods and technologies.

The innovative features and the outstanding optical quality provide the best image for the user.

The new generation of fibre-optic technology provides brilliant detail definition at the edges of the image as well as in the centre.

Recognition of even the smallest structures is thereby greatly facilitated. The natural colour rendition assists the user with better diagnostic.

The new generation is also designed for optimum hygiene.

Technical Data:**XQ 510/30 nopa nasopharyngoscope**

Field of view:	85°
Depth of field:	2,5 - 5 mm
Outer diameter insertion tube:	3,6 mm
Outer diameter distal end:	3,6 mm
Working length:	300 mm

Bending range of distal tip

Upward:	125°
Down:	125°