

Kaltlichtquelle
Cold Light Sources
Fuente de Luz Fria
Source à lumière froide
Sorgente Luce Fredda

**nopa® light XENON
Future 200**

180 Watt



XT 205/00 = 110-230 VAC/50/60 Hz

XT 204/00 Ersatzlampe
Spare bulb
Bombilla de recambio
Lampe de rechange
Lampadina di ricambio

Beschreibung und technische Daten auf der Rückseite
Description and technical data on reverse side

**nopa® light XENON 180 Watt****Beschreibung XT 205/00**

Die nopalight Xenon 180 Watt Lichtquelle kann bei allen Einsätzen angewendet werden, besonders dann, wenn eine starke Ausleuchtung oder eine Videoaufzeichnung notwendig ist.

Alle modernen Autoshutter-Videocameras werden mit ausreichender Lichtzufuhr versorgt. Die hohe Lichtausgangsleistung sorgt für absolut weißes Licht >5600°K und somit für detailgetreue Farbwiedergabe.

Unsere spezielle Lampenkonstruktion mit geringster Wärmeentwicklung führt zu absolut kaltem Licht und verhindert eine thermische Beschädigung der Lichtleiter.

Weitere Vorteile: Geringe Luftverwirbelungen, geringe Geräuschentwicklung, umweltfreundliche Lampe, einfacher Lampenwechsel.

Technische Daten:

Elektr. Versorgung: 110-230 VAC/50/60 Hz
Weltnetzteil

Farbtemperatur: > 5600° K

Maße
(B x H x T): 355 x 110 x 320 mm

Gewicht: ca. 9,5 kg

Zubehör:

Lampe XT 204/00: 180 Watt Xenon

nopa® light XENON 180 Watt**Description: XT 205/00**

The nopalight Xenon 180 watt light source can be used for all occasions particularly where strong illumination or a video recording is necessary.

All modern video camera automatic shutters are supplied with sufficient light. The high light output ensures absolute white light >5600° K and thus detail fidelity in colour reproduction.

Our special lamp construction which has the lowest of heat generation leads to absolute cold light and prevents thermal damage to the light conductor.

Further advantages are low air turbulence, low noise level, environmentally friendly lamp and very easy bulb replacement.

Technical Data:

Power Supply: 110-230 VAC/50/60 Hz
world power supply

Colour Temperature: > 5600° K

Measurements
(W x H x D): 355 x 110 x 320 mm

Weight: approx. 9,5 kg

Accessories:

Bulb XT 204/00: 180 Watts Xenon