



nopa Cam Digital Camera System

Beschreibung XT 302/00 / XT 302/01

Das **nopa' Cam** Camera System ist eine wirtschaftliche Office-Kamera mit ausgezeichneten Eigenschaften. Beste Bildqualität, scharfe Bildzeichnung und Klarheit werden durch den 1/4" CCD-Chip gewährleistet.

Die optimale Bilddarstellung aufgrund eines hochempfindlichen Mikro-Linsen-Chip sowie einer äußerst schnellen Blende wird auch unter schwierigen Beleuchtungsbedingungen sichergestellt. Der Kamerakopf ist leicht, ergonomisch und einlegbar.

Die Kamera ist mit digitaler Kontrast und Farbsteuerung ausgestattet. Der automatische Weißabgleich wird durch Schalterdruck aktiviert und ist im Memory gespeichert.

Aufgrund des hohen Automatisierungsstandards der **nopa' Cam S** ist das Kamerasystem sehr bedienerfreundlich.

Technische Daten:

XT 302/00 PAL - XT 302/01 NTSC

CCD Sensor:	1/4" Chip
Video-Auflösung:	über 500 Linien
Minimale Beleuchtung:	0,9 Lux
Video Ausgang:	FBAS oder Y/C (S-VHS)
Elektr. Versorgung:	110/115 VAC - und 230/240VAC
Gewicht:	1 kg

Weiteres Zubehör:

XT 302/02 Kamera mit C-Mount Gewinde PAL
XT 302/03 Kamera mit C-Mount Gewinde NTSC
XT 306/00 Bildteiler

XT 361/40 Zoom Adapter 15 - 25mm
XT 361/50 Zoom Adapter 25 - 40mm
verwendbar mit XT 302/02 und XT 302/03

Description XT 302/00 / XT 302/01

The **nopa' Cam** Camera System is an economic office camera with excellent characteristics. The best possible picture quality, picture delineation and clarity are ensured by the 1/4" CCD chip.

Optimal picture presentation is ensured by a highly sensitive micro-lens chip with an extremely rapid diaphragm even under difficult lighting conditions. The camera head is light, ergonomic and immersible. Our camera is equipped with digital contrast and color control. The automatic white compensation is activated by press button and is stored in the memory.

The **nopa' Cam S** high level of automatic operation makes this camera system extremely easy to operate.

Technical Data:

XT 302/00 PAL - XT 302/01 NTSC

CCD Sensor:	1/4" chip
Video Resolution:	more than 500 lines
Minimal Illumination:	0,9 Lux
Video Output:	FBAS or Y/C (S-VHS)
Power Supply:	110/115 VAC - and 230/240VAC
Weight:	1 kg

Optional Accessories:

XT 302/02 camera with C-mount adapter PAL
XT 302/03 camera with C-mount adapter NTSC
XT 306/00 beam splitter

XT 361/40 Zoom adapter 15 - 25mm
XT 361/50 Zoom adapter 25 - 40mm
to use with XT 302/02 and XT 302/03