

MONITOR



XT 403/21

Beschreibung und technische Daten auf der Rückseite
Description and technical data on reverse side



MONITOR

MONITOR

Beschreibung XT 403/21

Der nopa Monitor in 21" wurde entwickelt für den medizinischen Dauereinsatz. Er ermöglicht die Betrachtung von Untersuchungen und Operationen und garantiert eine effiziente Diagnoseerstellung und Überwachung während endoskopischer Eingriffe.

Der Dreifachstandard PAL, SECAM, NTSC ist automatisch vorwählbar.

Die hohe horizontale Auflösung von über 700 Linien gewährleisten optimale Bildqualität mit maximaler Detailtreue in Kombination mit den nopa-Kamerasystemen.

Eingebauter Lautsprecher und On-screen-Menü bieten einfache und vielfältige Anwendung.

Technische Daten:

XT 403/21	Farb-Monitor 52 cm/21" 100 - 240VAC, 50 - 60 Hz
Farbsystem:	PAL, SECAM, NTSC automatische Einstellung
Bildröhre:	Black matrix tube
Farbauflösung:	über 700 Linien horizontal
Video-/Audio-Eingang:	2 Composite Video-/ Audio-Eingänge Y/C-Eingang 1 VTR Video-/ Audio-Eingang RGB-Sync. IN/OUT Eingang

Maße: 21"

Höhe: 390 mm
Breite: 480 mm
Tiefe: 512 mm

Gewicht: 25,5 kg

MONITOR

Description XT 403/21

The nopa 21" monitor is developed for long period continuous use in medicine. It enables examinations and operations to be observed and guarantees an efficient diagnosis and surveillance during operations performed with an endoscope.

The triple standard PAL, SECAM and NTSC can be pre-selected automatically.

The high horizontal definition from more than 700 lines guarantees an optimal picture quality with maximum detail fidelity in combination with the nopa camera systems.

Integral loudspeakers and an On-Screen menu offer a simple and varied use.

Technical Data:

XT 403/21	colour monitor 52 cm/21" 100 - 240VAC, 50 - 60 Hz
Colour system:	PAL, SECAM, NTSC automatic selection
Screen:	Black matrix tube
Resolution:	more than 700 lines horizontal
Video-/Audio-Input:	2 Composite Video-/ Audio-Inputs Y/C-Input 1 VTR Video-/ Audio-Input RGB-Sync. IN/OUT

Measurements: 21"

Height: 390 mm
Width: 480 mm
Depth: 512 mm

Weight: 25,5 kg