

Flachbildschirm mit Standfuß
Flatscreen with base



XT 404/19

Beschreibung und technische Daten auf der Rückseite
Description and technical data on reverse side

**Flachbildschirm mit Standfuß**
Flatscreen with base**MONITOR****Beschreibung XT 404/19**

Der nopa Flachbildschirm in 19" wurde entwickelt für den Einsatz im Bereich der minimalinvasiven Chirurgie. Aufgrund der Aluminium Leichtbauweise können die Monitore an alle OP Arme, Trolleys und OP Lampen montiert werden.

Ein dichtes Gehäuse und eine Folientastatur bilden die Voraussetzungen, dass die Monitore feucht gereinigt werden können.

Die hohe Auflösung von 1280 x 1024 Pixel gewährleistet optimale Bildqualität mit maximaler Detailtreue in Kombination mit den nopa Kamerasystemen.

Der Standfuß ist höhenverstellbar und hat ein Kugelgelenk wodurch der optimale Blickwinkel erreicht werden kann.

Technische Daten:

XT 404/19	Medizinischer LCD Monitor 19" 100 - 240VAC, 50 - 60 Hz
Farbsystem:	PAL, SECAM, NTSC automatische Einstellung
Blickwinkel:	170° / 170°
Auflösung:	1280Px. x 1024Px.
Video-/Audio-Eingang:	DVI-I, VGA (über DVI), FBAS, S-VHS

Maße: 19"

Höhe: 386 mm
Breite: 442 mm
Tiefe: 67 mm

Gewicht: 6,9 kg

Optional:

XT 404/01 Standfuß allein
XT 404/02 Haltearm für nopa mobil

MONITOR**Description XT 404/19**

The nopa 19" flat screen was developed for use in the area of minimal invasive surgery. Due to the aluminium lightweight construction, the monitors can be mounted to all OP arms, trolleys and OP lamps.

A leak-tight housing and foil keyboard are required for wet-cleaning the monitors.

The resolution of 1280 x 1024 pixels guarantees optimum image quality with maximum detail precision in combination with the nopa camera systems.

The monitor stand can be adjusted (height) and has a ball joint, which allows the viewing angle to be optimally reached.

Technical Data:

XT 404/19	Medical LCD monitor 19" 100 - 240VAC, 50 - 60 Hz
Colour system:	PAL, SECAM, NTSC automatic selection
Field of view:	170° / 170°
Resolution:	1280Px. x 1024Px.
Video-/Audio-Input:	DVI-I, VGA (über DVI), FBAS, S-VHS

Measurements: 19"

Height: 386 mm
Width: 442 mm
Depth: 67 mm

Weight: 6,9 kg

Optional:

XT 404/01 Base only
XT 404/02 Swinging arm for nopa mobil