

Flachbildschirm mit Armsystem
Flatscreen with swinging arm system



XT 405/19

Beschreibung und technische Daten auf der Rückseite
Description and technical data on reverse side

**Flachbildschirm mit Armsystem**
Flatscreen with swinging arm system**MONITOR****Beschreibung XT 405/19**

Der nopa Flachbildschirm in 19" wurde entwickelt für den Einsatz im Bereich der minimalinvasiven Chirurgie. Aufgrund der Aluminium Leichtbauweise können die Monitore an alle OP Arme, Trolleys und OP Lampen montiert werden.

Ein dichtes Gehäuse und eine Folientastatur bilden die Voraussetzungen, dass die Monitore feucht gereinigt werden können.

Die hohe Auflösung von 1280 x 1024 Pixel gewährleistet optimale Bildqualität mit maximaler Detailtreue in Kombination mit den nopa Kamerasystemen.

Das Armsystem kann in beliebiger Höhe am nopa mobil montiert werden und bietet somit die bestmögliche Einstellung des Blickwinkels.

Technische Daten:

XT 405/19	Medizinischer LCD Monitor 19" 100 - 240VAC,50 - 60 Hz
Farbsystem:	PAL, SECAM, NTSC automatische Einstellung
Blickwinkel:	170° / 170°
Auflösung:	1280Px. x 1024Px.
Video-/Audio-Eingang:	DVI-I,VGA (über DVI), FBAS, S-VHS

Maße: 19"

Höhe: 386 mm
Breite: 442 mm
Tiefe: 67 mm

Gewicht: 6,9 kg

MONITOR**Description XT 405/19**

The nopa 19" flat screen was developed for use in the area of minimal invasive surgery. Due to the aluminium lightweight construction, the monitors can be mounted to all OP arms, trolleys, and OP lamps.

A leak-tight housing and foil keyboard are required for wet-cleaning the monitors.

The resolution of 1280 x 1024 pixels guarantees optimum image quality with maximum detail precision in combination with the nopa camera systems.

The arm system can be mounted at any height on the nopa mobil, and allows the best possible viewing angle to be set.

Technical Data:

XT 405/19	Medical LCD monitor 19" 100 - 240VAC,50 - 60 Hz
Colour system:	PAL, SECAM, NTSC automatic selection
Field of view:	170° / 170°
Resolution:	1280Px. x 1024Px.
Video-/Audio-Input:	DVI-I,VGA (über DVI), FBAS, S-VHS

Measurements: 19"

Height: 386 mm
Width: 442 mm
Depth: 67 mm

Weight: 6,9 kg