

VIDEO-RECORDER



XT 424/00 PAL
XT 424/01 NTSC

Beschreibung der technischen Daten auf der Rückseite
Description and technical data on the reverse side

Abbildungen und technische Änderungen vorbehalten
Change of model or technical data reserved



VIDEO-RECORDER

**VIDEO-RECORDER XT 424/00 PAL
XT 424/01 NTSC****Beschreibung: XT 424/00 - XT 424/01**

Der Videorecorder wurde speziell für den medizinischen Einsatz entwickelt. Die Wiedergabe erfolgt im S-VHS Format und ermöglicht eine hochwertige Aufnahmequalität sowie eine detailgetreue Bildwiedergabe. Ein integriertes Jog- Shuttle-Suchrad ermöglicht die einfache Lokalisierung gewünschter Bilder.

Neben einer großen Auswahl an Wiedergabe- und Suchlauffunktionen bietet das Gerät weitere Vorteile, z. B. menügeführte Funktionseinstellung, automatische Video-Kopfreinigung und Fehlercodeanzeige. Der Einsatz spezieller Videoköpfe sichert hochwertige und rauscharme Signale und einen guten Detailkontrast. Durch das kompakte und ergonomische Design ist das Gerät sehr benutzerfreundlich.

Eine Fernbedienung ist optional anschließbar. Das Gerät ist MedGV geprüft.

Technische Daten:

XT 424/00	S-VHS Video-Recorder, PAL 120 - 240VAC, 50 - 60 Hz
XT 424/01	S-VHS Video-Recorder, NTSC 120 - 240VAC, 50 - 60 Hz
Farbauflösung:	über 400 Linien horizontal
Videoein-/Ausgänge:	BNC 1 Vss 75 Ohm sync. neg. S-Video: 4-pin DIN Y: 1 Vss 75 Ohm sync. negativ C: 0,3 Vss Burstsinal, 75 Ohm sync. neg.
Audioeingänge:	Line In: CH1/CH2 Mikrofon CH2
Audioausgänge:	Line Out: CH1/CH2 Monitor Out Kopfhörerbuchse
Maße:	Höhe: 131,5mm Breite: 270mm Tiefe: 365,5mm Gewicht: 7,0kg

Weiteres Zubehör:

XT 424/10 Fernbedienung

Weiteres Zubehör auf Anfrage.
Technische Änderungen vorbehalten.

**VIDEO-RECORDER XT 424/00 PAL
XT 424/01 NTSC****Description: XT 424/00 - XT 424/01**

The video recorder has been specially developed for medical use. Reproduction is in S-VHS format and enables high-grade recording quality as well as detail fidelity picture reproduction. An integrated Jog and Shuttle search wheel enables desired pictures to be located easily.

Apart from a wide selection of play and search functions, the unit offers further advantages, such as menu guided function setting, automatic video head cleaning and an error code display. Use of special video heads ensure high grade and low noise signals and good contrast for details. The unit's compact and ergonomic design make it very easy to operate.

A remote control unit is available as an optional extra. The unit is approved according to the MedGV (German regulations for medical devices).

Technical Data:

XT 424/00	S-VHS Video-Recorder PAL 120 - 240VAC, 50 - 60 Hz
XT 424/01	S-VHS Video-Recorder, NTSC 120 - 240VAC, 50 - 60 Hz
Resolution:	more than 400 lines horizontal
Video-In-/Output:	BNC 1 Vss 75 Ohm sync. neg. S-Video: 4-pin DIN Y: 1 Vss 75 Ohm sync. neg. C: 0,3 Vss burst-level, 75 Ohm sync. neg.
Audio-Input:	Line In: CH1/CH2 Mikrophone CH2
Audio-Output:	Line Out: CH1/CH2 Monitor Out Headphone
Measurements:	Height: 131,5mm Width: 270mm Depth: 365,5mm Weight: 7,0kg

Optional Accessories:

XT 424/10 Remote Control

Further Accessories upon request
Technical specifications subject to change without notice.