

VIDEO-PRINTER



XT 440/00 PAL
XT 440/01 NTSC

Beschreibung der technischen Daten auf der Rückseite
Description and technical data on the reverse side



VIDEO-PRINTER

VIDEO-PRINTER

XT 440/00 PAL
XT 440/01 NTSC

Beschreibung: XT 440/00 - XT 440/01

Der Sony Videoprinter UP 2100P arbeitet im Farb-thermosublimations-Verfahren und ist sehr anwenderfreundlich. Fotorealistische Ausdrücke werden aufgrund hochauflösender Thermoköpfe sowie extrem schneller Druckwalzentechnik realisiert.

Eine Fernbedienung ist optional erhältlich.

Technische Daten:

XT 440/00	Video-Printer, PAL 220 - 240VAC, 50 - 60 Hz
XT 440/01	Video-Printer, NTSC 220 - 240VAC, 50 - 60 Hz
Thermokopf: Farbabstufung:	608 Elemente, 6,72 Punkte/mm 24 bit, 256 Abstufungen über 16,7 Mio Farben
Bildgröße:	80 x 110mm
Bildpunkte:	588 x 720 dpi
Videoein-/Ausgänge:	FBAS (BNC) CCIR 1 Vp-p 75 Ohm sync. negativ Y/C 4-pin Mini DIN Y: 1,0Vp-p 75 Ohm sync. negativ C: 0,3Vp-p 75 Ohm sync. negativ

Maße:	Höhe:	125mm
	Breite:	370mm
	Tiefe:	417mm
	Gewicht:	9,0kg

Weiteres Zubehör:

XT 440/10	Fernbedienung
XT 440/50	Papier und Farbband

Weiteres Zubehör auf Anfrage.
Technische Änderungen vorbehalten.

VIDEO-PRINTER

XT 440/00 PAL
XT 440/01 NTSC

Description: XT 440/00 - XT 440/01

The Sony video printer UP 2100P employs the colour thermo-sublimation process and is very simple to use. Because of the high definition thermo-heads and extremely fast pressure roller technology, prints in photographic quality are produced.

A remote control unit is available as an optional extra.

Technical Data:

XT 440/00	Video-Printer, PAL 220 - 240VAC, 50 - 60 Hz
XT 440/01	Video-Printer, NTSC 220 - 240VAC, 50 - 60 Hz
Thermal head: Colour graduation:	608 elements, 6,72 dts./mm 24 bit, 256 graduations more than 16,7 Mio colours
Picture size:	80 x 110mm
Pixel:	588 x 720 dpi
Video-In-/Output:	FBAS (BNC) CCIR 1 Vp-p 75 Ohm sync. neg. Y/C 4-pin Mini DIN Y: 1,0Vp-p 75 Ohm sync. neg. C: 0,3Vp-p 75 Ohm sync. neg.

Measurements:	Height:	125mm
	Width:	370mm
	Depth:	417mm
	Weight:	9,0kg

Optional Accessories:

XT 440/10	Remote Control
XT 440/50	Paper and color ribbon

Further Accessories upon request
Technical specifications subject to change without notice.