



XT 448/00

Beschreibung und technische Daten auf der Rückseite
Description and technical data on reverse side

**MEDICAPTURE****Beschreibung XT 448/00**

MEDICapture ist ein neuer Weg um Videobilder auf einem USB - Datenträger zu speichern. Verschiedene Videosysteme und Eingänge machen das MEDICapture zu einem sehr flexiblen Instrument für die Inspektion und Dokumentation.

MEDICapture kann an jeder Kamera, jedem Videoendoskop und Mikroskop mit analogem Ausgang betrieben werden. Bilder werden in Standardformaten (JPEG/BMP) mit automatisch generierten Dateinamen abgespeichert.

Das MEDICapture generiert ebenfalls eine Vorschauansicht zur Überprüfung oder Löschen der Dateien auf dem Monitor. Die Dateien können sehr einfach auf jeden Computer mit USB Anschluss übertragen werden.

Technische Daten:

XT 448/00	MEDICapture 100 - 240VAC, 50 - 60 Hz
Farbsystem:	PAL, NTSC einstellbar
Videoeingang:	1 Composite Videoeingang 1 S-VHS / Y/C-Eingang
Videoausgang:	1 Composite Videoausgang 1 S-VHS / Y/C-Ausgang
Medium:	USB - Speichermedium
Dateiformat:	DICOM, DICOM Hochauflösend JPEG, BMP
Bildgröße:	800 x 600 Bildpunkte
Standard:	FCC, UL
Höhe:	63 mm
Breite:	240 mm
Tiefe:	200 mm
Gewicht:	680 g

Weiteres Zubehör:

XT 448/10	Fußschalter
-----------	-------------

MEDICAPTURE**Description XT 448/00**

MEDICapture is the new simple way to capture video images on a USB - Memory-Stick. Different video systems and different input standards make the MEDICapture a very flexible instrument for inspection and documentation.

MEDICapture can be connected to any camera, video-scope and microscope with analog outputs. Images are compressed and stored with auto generated filenames and information in standard file formats (JPEG/BMP).

MEDICapture generates a preview on the display to verify, erase or show the images. The files can be uploaded very simple to every computer with USB connection.

Technical Data:

XT 448/00	MEDICapture 100 - 240VAC, 50 - 60 Hz
Colour system:	PAL, NTSC switchable
Video input:	1 Composite Video input 1 S-VHS / C-Input
Video output:	1 Composite Video output 1 S-VHS / C-output
Media:	USB - Stick
File Format:	DICOM, DICOM True color JPEG, BMP
Image size:	800 x 600 Pixel
Approvals:	FCC, UL
Heigh:	63 mm
Width:	240 mm
Depth:	200 mm
Weight:	680 g

Optional Accessories:

XT 448/10	Foot-Switch
-----------	-------------