

CD Video documentation

nopa® Cap

CD Video documentation



PAL / NTSC

XT 450/00

=

XT 310/00 Netzkabel / mains cord

XT 450/02 Fernsteuerungskabel / remote control cord

XT 450/03 CD-RW Rohling / CD-RW recordable

XT 450/04 Software Pack

XT 344/02 Kabel / cable S-VHS

XT 346/00 Kabel / cable BNC

Beschreibung und technische Daten auf der Rückseite

Description and technical data on reverse side



CD Video documentation

Beschreibung XT 450/00

nopa[®] Cap motion ist die neueste Video- und Bildspeichertechnik auf dem Markt. Die intraoperative "Ein-Knopf" Bedienung garantiert den unkomplizierten Gebrauch für den Chirurgen oder das OP-Team, ohne dass zusätzliche Computerkenntnisse oder Personal benötigt wird.

Günstige, haltbare und universell verwendbare Speichermedien (CD-R und CD-RW) ermöglichen raumsparende Archivierung von Daten unabhängig von bereits vorhanden In-Haus-Standards.

Das Editieren, Speichern und die Verbindung aller Bilder und Videodaten mit anderen Applikationen ist mit jedem gängigen PC außerhalb des OP's möglich. Dadurch wird kein zusätzliches Personal im OP gebunden.

Technische Daten:

XT 450/00

Viedo-Standards:	PAL und NTSC
CD-Sampling-Rate:	44.1 kHz 7 16 Bit
Video-Eingänge:	1xY/C (S-VHS) 1x Composite Video (BNC)
Video-Ausgänge:	1xY/C (S-VHS) 1x Composite Video (BNC)
Elektrische Versorgung:	100/115 VAC und 230/240 VAC
Gewicht:	6,0 kg

Weiteres Zubehör:

XT 450/10	Fußschalter
XT 450/11	Seriellles Verbindungskabel

Description XT 450/00

nopa[®] Cap motion is the latest technics in digital video- and image-acquisition and storage. The intra-operative single-button-control guarantees unprecedented ease of use for the surgeon without complicated selection of formats or patient data. No special computer knowledge or additional OR staff are needed to control this device.

Cost effective, durable and universal standards storage media (CD-R and CD-RW) allow for space-saving recording, independent from other in-house technical standards.

Editing, storing and integrating the image and video data into other applications is possible on any PC outside the OR without using valuable space and staff in the operating room.

Technische Data:

XT 450/00

Viedo-standards:	PAL and NTSC
CD sampling rate:	44.1 kHz 7 16 Bit
Video input:	1xY/C (S-VHS) 1x composite video (BNC)
Video output:	1xY/C (S-VHS) 1x composite video (BNC)
Power supply:	100/115 VAC and 230/240 VAC
Weight:	6,0 kg

Optional accessories:

XT 450/10	foot pedal for remote control function
XT 450/11	direct link cable