


Suction Irrigation System

nopa®

Suction Irrigation System



XU 310/00

Beschreibung und technische Daten auf der Rückseite
Description and technical data on reverse side



Suction Irrigation System

nopa® Saug/Spül System

Beschreibung XU 310/00

- Höchstes Leistungsniveau
- Absaugen/Spülen mit bis zu 3 l/min
- Großer Absaugbehälter
- Standby-Funktion für minimalen Verschleiß
- Optimalste Sterilbedingungen durch spezielle Filter
- Einfaches, ergonomisches Bedienkonzept
- Variable Flaschenposition (am Gerät oder Gerätewagen)
- Keine Einstellarbeiten – Saug- und Spüleistung individuell am Instrument regelbar
- Integrierter Überlaufschutz
- Geringer Geräuschpegel durch Stand-by-Funktion
- Qualitativ hochwertige Bauteile

Technische Daten:

Netzspannung:	100 – 115 Vac 220 – 240 Vac
Leistungsaufnahme:	120 VA
Frequenz:	50 / 60 Hz
Spüleistung:	bis 3000 ml/min
Schutzklasse:	I
Nach Med.-Produkte-Richtl.:	CE 0123
Gewicht:	ca. 4,7 kg
Abmessung: H x B x T	14,5 x 30,0 x 32,0 cm

XU 310/00

Saug/Spül System inklusiv:

XU 310/91	Schlauchset mit LL-Anschlüssen
XU 310/92	Einstechrohr
XU 310/93	2 x Flaschenhalter
XU 310/94	Mikrofilter
XU 310/95	Netzkabel
XU 310/96	2 Liter-Absaug-Behälter inkl. Saugdeckel

Technische Änderungen vorbehalten

nopa® Suction Irrigation System

Description XU 310/10

- Highest efficiency
- Aspirates and irrigates at up to 3 l/min
- Large two-liter aspiration container
- Standby function minimizes wear and tear
- Special filters ensure sterility
- Simple ergonomic operation
- Flask may be attached to unit or carrier cart
- No special adjustment necessary. Aspiration and irrigation are regulated and selected by the instrument controls
- Integrated overflow protection
- Stand by function for low noise level
- High-quality components

Technical Data:

Line voltage:	100 – 115 Vac 220 – 240 Vac
Power consumption:	120 VA
Frequency:	50 / 60 Hz
Rinsing performance:	up to 3000 ml/min
Protection class:	I
Acc. To MDD:	CE 0123
Weight:	approx. 4,7 kg
Dimension: H x W x D	14,5 x 30,0 x 32,0 cm

XU 310/00

Suction Irrigation System including:

XU 310/91	tubing set with LL connections
XU 310/92	puncture tube
XU 310/93	2 x bottle holders
XU 310/94	microfilter
XU 310/95	power cord
XU 310/96	2 liter container incl. suction lid

Technical specifications subjects to change without notice.